

西安瑞禧生物科技有限公司 Xi'an city ,Shaanxi Province,China Web:http://www.xarxbio.com Tel: 0086 29 86354885 E-mail:sales@xarxbio.com

Certificate of Analysis

Product name:	Concanavalin A-C	y5
---------------	------------------	----

Product code: CN5Y

Product lot: 2503CN5Y3

Source: purified from Canavalia ensiformis

Maximum absorption: 650 nm

Maximum emission: 670 nm

Modification level: 0.8 Cy5 per con A monomer

Concentration: 0.4 mg/ml

Preservative: none

Solvent: water

Storage conditions: -Aliquot and store at -20 ℃ or colder

- -Concentrated concanavalin A agreggates quickly at room temp and slowly at 4°C.
- -Freezing and thawing should be done quickly to avoid high localized concentrations which facilitate aggregation and precipitation.

This product is for research use only, Not for human or animal use