

Certificate of Analysis

Product name Product code Product lot	Streptavidin Thiol SAVT 2602SAVT9
Source:	Recombinant protein produced in E.coli
Protein concentration:	1 mg/ml
Thiols/monomer:	1.5 thiols per monomer
Additives:	20 mM sodium citrate pH 4.00 3 mM EDTA
Thawing:	Do NOT thaw in refrigerator, since this is too slow and may create excessively high local protein concentration in the portion that thaws first. Thaw at room temperature or speed the process by squeezing in your hand. Make sure to return to 4 C immediately after thawing.
Storage conditions:	KEEP COLD. Store at 4 C for up to 1 week Freeze at -20 or colder for long term storage. Avoid multiple freeze/thaw cycles.

This product is for research use only, Not for human or animal use